

STATSEAL[®]

StatSeal[®] is intended for use under the care of a healthcare professional for the temporary external control of bleeding from wounds and procedures.

Each application includes one (1) bottle of **StatSeal Topical Powder**, one (1) powder containment device (PCD), one (1) PCD foam insert and one (1) product description label.

Instructions for Use

1. Open sterile packaging and remove contents. Remove paper backing from PCD to expose adhesive.
2. Center PCD around site with adhesive side down.
3. Fill PCD to brim with powder; **do not overfill***. Ensure powder completely covers and surrounds site. Place PCD foam insert over powder to cover.
4. Hold **continuous** pressure over site for 2 minutes. Extend hold times as necessary for anticoagulated patients.
5. Release pressure slowly by rolling fingers off of PCD.
6. Place transparent dressing per protocol. Apply product description label over dressing.
7. Leave powder in place until next scheduled dressing.

*The bottle contains extra powder. If bleeding continues, reapply more powder as necessary with manual compression until bleeding stops.

Additional instructions and procedural tips are available at:

www.statseal.com

INS-SS-04

STATSEAL®

Ingredients: Hydrophilic polymer and potassium ferrate

WARNINGS: FOR EXTERNAL USE ONLY. Not intended to replace stitches. Do not apply to visibly infected wounds. May be harmful if inhaled or ingested. Avoid aspirating product into the nasal or oral cavity. Avoid contact with eyes; if contact should occur, thoroughly rinse eyes with water. Keep out of reach of children.

CAUTION: Store in a clean, dry place. Use by expiration date listed. Discard if product is damaged, opened or if powder is clumped and does not flow freely. Discard contents after use; powder is guaranteed to work only if it has NOT been left exposed to air and humidity. After use, the remaining bottle contents and the applicator tip are no longer sterile.

**The symbols glossary is provided electronically at:
www.biolife.com/symbols**



STERILE R



Biolife, LLC

8163 25th Court East, Sarasota, FL 34243, USA

+1 941.360.1300 • customer.care@biolife.com



0543



EMERGO EUROPE

Prinsessegracht 20, 2514 AP The Hague, Netherlands

Patents 6,521,265, 6,187,347 and 8,007,833.

Other U.S. and international patents pending.